

IHE Integration Statement

Vendor	Product Name	Version	Date
Orion Health Inc.	Clinical Portal	6.7	January 25, 2011

This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:

Integration Profiles Implemented	Actors Implemented	Options Implemented
Patient Identifier Cross Referencing for MPI (PIX)	Patient Identifier Cross-Reference Consumer	PIX Update Notification
Patient Identifier Cross Reference HL7 V3 (PIXv3)	Patient Identifier Cross-Reference Consumer	PIXv3 Update Notification Transaction
Patient Demographics Query (PDQ)	Patient Demographics Consumer	Patient Demographics and Visit Query
Patient Demographics Query HL7 V3 (PDQV3)	Patient Demographics Consumer	None
Consistent Time (CT)	Time Client	None
Audit Trail and Node Authentication (ATNA)	Secure Node	None
Cross-Enterprise Document Reliable Interchange (XDR)	Document Source	None
Cross-Enterprise Document Sharing (XDS.b)	Document Consumer	None
Cross-Enterprise Document Sharing (XDS.b)	Document Source	None
Cross-Enterprise Sharing of Scanned Documents (XDS-SD)	Content Consumer	View Option
Cross-Enterprise Sharing of Medical Summaries (XDS-MS)	Content Consumer	View Option
Exchange of Personal Health Record Content (XPHR)	Content Consumer	View Option



IHE Integration Statement

Vendor	Product Name	Version	Date
Orion Health Inc.	Clinical Portal	6.7	January 25, 2011

This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:

Healthcare Information Technology Standards Panel (HITSP)

HITSP Consistent Time (T16)	Time Client	None
HITSP Collect and Communicate Audit Trail (T15)	Secure Node	None
HITSP Secured Communication Channel Transaction (T17)	Secure Node	None
HITSP Manage Sharing of Documents (XDS.b Option)	Document Consumer	None

Internet address for vendor's IHE information

<http://www.orionhealth.com/concerto/ihe>

Links to Standards Conformance Statements for the Implementation

HL7
DICOM

Links to general information in IHE

In North America: <http://www.rsna.org/IHE>
 In Europe: <http://ihe-europe.org>
 In Japan: <http://jira-net.or.jp/ihe-j>

